Cheatography

Biocompatibility of Dental Materials Cheat Sheet by Carm (Carmilaa) via cheatography.com/49544/cs/15294/

Definition:

Ability of a material to elicit an appropriate biological response in a given application in the body.

Placement of a material creates an interface that is not normally present.

Interface = Not static, includes dynamic interactions

Activity of interface depends on: location of material, duration in body, properties of material and health of host.

Key concepts of Biocompatibilty:

Biomaterials are not biologically inhert:

There are no inhert materials. Interactions depend on the material, the host and forces/conditions placed on the material. Material affects the host and host affects the material.

Dynamic Process:

Body may change due to aging/disease. Material may change due to corrosion/fatigue. Changes may alter conditions where initially an appropriate/desired biological response took place. Biological response to a material is an ongoing process.

Property of a Material and Its Environment:

Interaction

Organisations That Test Dental Materials:

> Food and Drug Administration
 (FDA)
 > American National Standards
 Institute

> American Dental Association
(ADA)

 > International Organization for Standardation

Measuring Biocompatibility of Materials:

Location:

Important to it's overall biological response. Determines if material will be covered by soft or mineralized tissue, will it be external to oral epithelium or if it will communicate through epithelium and will the material be directly exposed to bone, tissues , blood and saliva or will there be a barrier (such as enamel or dentin).

Duration:

More stringent tests for longer duration.

Stress:

Stresses placed on material (physical, chemical, thermal or occlusal or reaction to salivary proteins) may react unfavourably.

In-vitro Testing:

Performed outside an organism

First screening test done in test tubes, cell culture dish, flask.

Material is placed in direct/indirect contact with some biological organism.

Advantages:

> Fast

> Inexpensive

> Easily standardised

> larger scale screening conditions can be tightly controlled to provide the highest quality of scientific

rigor.

Disadvantages:

> Potential lack of relevance to the in vivo use of material.

> Lacks the complex co-ordination of systems (immune system, inflammatory system, circulatory system)

Animal Testing:

Place a material into contact with an intact organism

Distinct difference, exposure without regard to the material's final use

Advantage: Ability to allow an intact biological system to respond to a material.

Disadvantages:

> interfere with the many complex biological systems

- > expensive
- > difficult to control
- > ethical considerations
- > questionable

Usage/Invitro Testing:

Material is placed in an environment, clinically relevant to the use of the material in clinical practice

Humans (clinical trial) and Animals closely resembling humans

Advantage: Gold standard

Disadvantage:

- > Expensive
- > Approval
- > Legal liabilities
- > Time consuming

Adverse Affects from Dental Materials:

1. Toxicity:

normally the first screening test.

2. Inflammation:

second fundamental type of biological response resulting either from allergy/ toxicity, first infiltration of neutrophils, then monocytes and other lymphocytes.

3. Allergic Reaction:

the most common to lay people. **Type I**- refers to an immediate atopic /anaphylatic reaction when an antigen reacts with mast cell/basophils. **Type II**- is a cytotoxic hypersensitivity reaction. **Type III**-immunecomplex hypersensitivity reaction. **Type IV** -delayed or cell-mediated hypersensitivity. **Type V** - stimulating-antibody reaction. **Type VI**- anti- body dependent, cell mediated cytoxicity reaction.

By C chear

By **Carm** (Carmilaa) cheatography.com/carmilaa/

Published 3rd April, 2018. Last updated 3rd April, 2018. Page 1 of 2.

Sponsored by ApolloPad.com

Everyone has a novel in them. Finish Yours! https://apollopad.com

Cheatography

Biocompatibility of Dental Materials Cheat Sheet by Carm (Carmilaa) via cheatography.com/49544/cs/15294/

Adverse Affects from Dental Materials: (cont)

4. Mutagenicity:

results when the components of a material alter the basepair sequences of DNA in cells.

Local and Systemic Effects:

Local: Pulp of the tooth, periodontium, root apex, soft tissues (buccal mucosa/tongue)

Systemic: Ingestion,

absorption in the gut, , inhaled vapor, release a the tooth apex, absorption through the o/mucosa. As well as, Systemic biological response depends on; duration and concentration of the exposure, excretion rate of the substance, site of exposure.

Principles of Adverse Effects:

Degradation process: The biological response to the corrosion products depends on the amount, composition, and form of these products as well as their location in the tissues.

Surface characteristics:

research has shown that for all dental materials , the surface is very different compared to the interior.

Examples:

Latex:(3 categories: irritant, allergic & type1 hypersensitivity reaction) In dentistry; rubber dam, gloves, bite blocks etc. Reactions to latex ; localized rashes, swelling to a more serious wheezing and anaphylaxis. Dermatitis of the hands (most common).

Principles of Adverse Effects: (cont)

Mercury&Amalgam: No data to show that mercury released from dental amalgam is harmful.

Estrogenicity: Ability of a

chemical to act as the hormone estrogen in the body The placement of composites has been questioned. Xenoestrogen bisphenol A (BPA) synthetic starting point for all Bis-GMA composites fear is that the release of these substances might alter normal cellular development or maintenance if the BPA has estrogenic effects.

Nickel: Common component of many dental alloys including those used for crowns, fixed partial dentures, removable partial dentures and some orthodontic appliances. Endodontic files Reactions are often subtle, resemble periodontal inflammation Occurs primarily outside the mouth

Beryllium:

if inhaled cause chronic inflammatory conditionberylliosis, the lung is engorged with lymphocytes and macrophages. T cells in susceptible individuals proliferate locally in the lung tissue ,in delayed hypersensitivity reaction too the beryllium metal. Occurs in individuals with a hypersensitivity to beryllium and may occur from inhalation of beryllium dusts, salts, fumes as those encountered when casting beryllium-containing alloys.

Adver	se Effects of Restorative Material:
TABLE 5:	Biocompatibility considerations of various dental restorative materials ^{60,11}
Restorative Material	Biocompatibility Consideration
Dental Amalgam:	 No adverse pulpal responses from mercury Correston may limit marginal backeys, but in the long term may load to breakdown of marginal integrity, especially with low copper amalgums Licktened resources reported Thermal conduction to pulp Mercury allergy (95)
Resin-Base Composites	Documented extragancing issue Very little reach on systemic bicompatibility Allergic to resin composite ingredients (8%) Incomplete polymerization leading to degadation, tracking, and imperfect bonding Predisposed to polymerization shrinkage Associated with adverse load pulguland dentin reactions, development of recorrect carles, and pain Higher proprietion of streptococcus mutans leading to secondary carles
Glass lonomer Cements	 Fee documented systemic adverse effects. Early pulger actions, although less than with cements or composite resins, and with rapid recovery. Hydraulic pressure and etchicing during placement may irritate the pulp. Goad affesion, minimal leskage at margins, high biocompatibility. Laaching of component materials dires soperunity for solve release of thoride
Gold Foil and Cast Alloys	Inert: sensitivities are rare Potential pulpal reactions due to condensation Gold contact allergy (23%)
Ceramics	Inert material No long-term data on biocompatibilility Possibility of silica granulomas

Sponsored by ApolloPad.com

Everyone has a novel in them. Finish Yours! https://apollopad.com



Published 3rd April, 2018. Last updated 3rd April, 2018. Page 2 of 2.