

Pregnancy / hcG

Controls

Brand Name	Expire/BUD	Storage
Quantamatrix UA controls	Printed date on box	36°F-45°F
SureVue	Printed date on box	36°F-45°F
Alere	BUD 30days (Rm temp) ; 90days (fridge)	36°F-45°F or Rm temp

Test Kits

Brand Name	QC frequency	Storage	Exp. Date
<i>Medline</i>	Lot number & every 30 days	56°F-86°	Date on box



By **mwaford**
cheatography.com/mwaford/

Not published yet.
Last updated 8th February, 2023.
Page 2 of 2.

Sponsored by **CrosswordCheats.com**
Learn to solve cryptic crosswords!
<http://crosswordcheats.com>

