



Lead Placement



Peripheral Pulse Locations



Peripheral Pulses Data

Peripheral Pulses

Rate: Number per Minute—60-90
Rhythm: Regular or Irregular
Symmetry: Even Uneven

Amplitude:
 4 = Bounding (fluid, thyroid, CRD, CHF)
 3 = Increased (MOST PT WISH)
 2 = Normal (Stable)
 1 = Weak (Dehydration, Occlusion)
 0 = Absent or No palpable (Acute)

If unable to assess pulse Doppler (US)
 Doppler use is a potential acute condition.

Capillary Refill: <3 Seconds Normal

C

By **dcox80**
cheatography.com/dcox80/

Not published yet.
 Last updated 21st October, 2024.
 Page 1 of 1.

Sponsored by **CrosswordCheats.com**
 Learn to solve cryptic crosswords!
<http://crosswordcheats.com>