

Cleaning

Add 2 tablespoon of cleaner to 15 ltr of water at 60C

Install malt tube and counterflow chiller

Circulate cleaning solution for 15 mins

Install recirculation tube

Recirculate cleaning solution for 5 mins

Scrub any remaining stains

Rinse thoroughly

Mash

Install malt pipe and heat to 67C Weigh out malt

Add malt to malt pipe

Start to do the Mash! Usually for 60 mins but may vary

After 40 mins of the Mash, start to heat the Sparge water

Enzyme Temperature Ranges

Enzyme	Temperature Range (C)
--------	-----------------------

Phytase	30-53
---------	-------

Beta-Glucanase	35-55
----------------	-------

Peptidase	45-53
-----------	-------

Proteinase	50-59
------------	-------

Beta-amylase	54-66
--------------	-------

Alpha-amylase	66-71
---------------	-------



By euuen
cheatography.com/euuen/

Not published yet.
Last updated 10th April, 2025.
Page 1 of 1.

Sponsored by **CrosswordCheats.com**
Learn to solve cryptic crosswords!
<http://crosswordcheats.com>