

### Lab 1

Accuracy: How close the measured value is to the "true value"

Precision: Fitness of a measurement.  
Reproducibility of a measurement.

Density: Mass/Volume

All sig figs should be used in measurements: i.e. 16.4400g not 16.44g

Molar Mass = Mass(g)/# of Moles

### Lab 3



By **Flxstr1**  
[cheatography.com/flxstr1/](https://cheatography.com/flxstr1/)

Not published yet.  
Last updated 28th November, 2022.  
Page 1 of 1.

Sponsored by **CrosswordCheats.com**  
Learn to solve cryptic crosswords!  
<http://crosswordcheats.com>